
Endodonzia *Endodontics*

Quality Products
Made in Germany

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Some of the products and designations mentioned in the text are trademarked, patented or copyrighted. The absence of a special reference or the sign ® should not be interpreted as the absence of legal protection.

This publication is copyrighted. All rights, also with regard to translation, reprint and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or processed using electronic systems in any form or by any means (photocopying, microfilm or other methods) without the written permission of the editor.

*Colours and products subject to alterations.
Printing errors excepted.*

As at September 2014

4 - 10	Apertura cavità d'accesso <i>Preparation of the access cavity</i>
11 - 19	Punte a vibrazione sonora <i>Sonic tips</i>
20 - 23	Preparazione manuale <i>Manual preparation</i>
24 - 25	Glide Path meccanico <i>Mechanical glide path</i>
26 - 33	Preparazione Meccanica <i>Mechanical preparation</i>
26 - 27	• F6 SkyTaper <i>F6 SkyTaper</i>
28 - 29	• F360® <i>F360®</i>
30 - 31	• EndoPilot mobil <i>EndoPilot mobil</i>
32 - 33	• E-Drive <i>E-Drive</i>
34 - 42	Otturazione <i>Obturation</i>
43 - 49	Accessori <i>Auxiliaries</i>

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Preparation of the access cavity

The preparation of the access cavity is the first and possibly the most important step of the entire endodontic treatment.

Having the right tools at hand, preferably in an efficient and well thought-out set, makes this task considerably safer and easier.

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H47L

Pear, long

Matches H47L carbide finisher series

		5	5	5	5	5	5	5	5	5
Misura - Size	Ø 1/10 mm	008	010	012	014	016	018	021	023	027
US No.		1SM	2SM	3SM	4SM	5SM	6SM	7SM	8SM	10SM

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
Round bur, made of ceramics

High-efficiency cutting design for excavating

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
Round

Staggered toothings with high-efficiency cutting design for excavating

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

For isthmus preparation, slim neck for improved vision

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

Combination bur with diamond coated working part and carbide tip, for lateral substance removal when preparing an access cavity

Apical problems solved the coronal way!

According to this motto for a successful endodontic treatment, the roof of the pulp chamber is completely removed. To create a straight access cavity to the root canal, lateral gain of space and the removal of overhanging material are particularly important. Both abrasive and cutting instruments are available for this purpose. They can be used individually or in combination, as required.

Instruments with safe, non-cutting ends: Diamond and tungsten carbide burs with safe ends are available to protect the base of the pulp chamber and the canal entrance.

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AlphaKite File Taper 10, length 15 mm or 19 mm

Universal reamer for the straight root canal entrance area, nickel-titanium with TiN coating

In countries other than Germany and Austria the packing unit is 5 instead of 6

new

LD 1250 A.000

Apertura Camera Pulpale Dr. Giuseppe Squeo-Bari
Opening the pulp chamber Dr. Giuseppe Squeo-Bari

●	S6830L.314.012	1	
	H269GK.314.012	1	
	H269QKG.314.016	1	
●	H1SM.314.014	1	
	H1.316.012	1	
●	OP10L19.204.030	1	

SonicLine

SonicLine

Orthograde preparation of the pulp chamber and preparation of the cervical third of the root canal as part of an endodontic treatment

Advantages:

- Quick preparation and removal of old root fillings
- Easier retrieval of root canals
- Enlargement of obliterated canals
- Preparation of straight canal access cavities without weakening the crown
- Useful for removing hard root fillings, cements or posts
- Controlled, gentle preparation without steps and protruding material
- Excellent vision

The sonic tips can be used in the air driven sonic hand piece SF1LM (air scaler) and reprocessed with the help of a rinse adapter in a Miele washer/disinfector.

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Large bud

Initial preparation of the access cavity and removal of protruding substance

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Tapered

Retrieval of fine and calcified canals, opening of the upper canal portions during revision

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Tapered

Alternative to the SF67 with a more pronounced angle

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Small bud

Finishing of the access cavity, minimal refreshing of the dentin layer and removal of residues of root canal fillings

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Tapered

Enlargement of long and wide canals, unblocking of fractured instruments, removal of root fillings made of gutta-percha and soft cements

SonicLine

Sonic tip for activating rinsing liquids during endodontic treatments.

Thanks to sound activated movements and micro currents, the efficiency of rinsing liquids inside the root canal is greatly increased, safely removing bacteria, residues of pulp tissue, loose dentin chips and the smear layer.

Advantages:

- Even more efficient rinsing of the root canal
- Made of highly flexible nickel titanium, with a titanium-nitride coating
- Small instrument diameter for rinsing narrow root canals
- The same instrument diameter for all canals, thus eliminating the need to change the instrument
- Safe end and absence of toothing to avoid inadvertent removal of substance from the canal wall
- Laser marks to indicate the depth

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




1

SonicLine

For retrograde endodontic treatments as part of an apicectomy

Clinical sequence:

Preparative surgical steps

1. Create the access through the jaw bone and work on the bone in the periapical region. Prepare a small bone window with a repositionable lid to ensure a proper preparation of the canal with the endodontic tips for retrograde treatments.

2. Remove the root apex in a right angle to the tooth axis, for example with the tungsten carbide bone cutter Komet® H254E

3. Remove any diseased tissue by means of a sharp spoon curette and stop the resulting bleeding

Retrograde root preparation with SonicLine endodontic tips for retrograde treatments

4. Retrieval of the root canal and exposure of the canal entrance with the sonic tips SF56 (curved to the left) or SF57 (curved to the right)

5. Preparation of the cavity up to a depth of 3 mm with the SF16 (curved to the left) or SF17 (curved to the right) In case of extremely curved root canals, it might be helpful to open and enlarge the canal with the specially shaped sonic tip SF55

6. Preparation of an undercut in order to prevent loss of the canal filling with the tip SF20 (curved to the left) or SF21 (curved to the right)

7. Finally, fill the root with a material suitable for retrograde fillings

Advantages:

- Minimally invasive treatment without need to prepare a large bone window
- Axial work, even in very crowded conditions
- Easier work thanks to double-angled tips
- Slender tips for excellent vision in all jaw regions
- Simplified preparation of undercuts for permanent retention of the retrograde root filling

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Instruments for root canal preparation, Hedstroem files

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

Instruments for root canal preparation, root canal

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
9107.654.S1

2 x 030
2 x 035
2 x 040

PathGlider for the mechanical preparation of a glide path

The preparation of a glide path entails creating smooth canal walls, along which the NiTi instrument can glide with ease in order to reach the working length and to prepare the apical foramen to a size that is compatible with that of the first rotary instrument to be used. This avoids the risk of a "taper lock", thus reducing the probability of the rotary instrument fracturing in the root canal.

The manual glide path is established with a steel K-file, which is one of the most delicate phases of the process. Inadvertent errors are most likely to occur during this phase, such as undesirable steps, via falsa and dentin plugs, which can lead to the failure of the entire endodontic treatment.

In order to eliminate this type of error, the glide path can nowadays be created with rotary NiTi instruments, provided that the working length has been reached with a K-file in size 010 beforehand.

The PG03 can establish a mechanical glide path in just one step!

The glide path significantly reduces the risk of rotary NiTi instruments fracturing in the root canal and is essential for getting to know the anatomy of the root canal.

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
PathGlider file, taper 03, length 25 mm, spiral shaped colour code, single-use file for the mechanical establishment of a glide path, sterile packed, nickel-titanium

F6 SkyTaper – A new level of flexibility

The first one-file system on the market that unites flexibility of the file with outstanding versatility of use.

- *Unlimited versatility. The F6 SkyTaper System is composed of highly flexible nickel-titanium files in five different sizes and three lengths, which offers the user the ideal solution for all types of root canal treatment.*
- *Thanks to their taper of .06, the F6 SkyTaper files are equally suitable for lateral, carrier based or thermoplastic obturations.*
- *Without restrictions. The F6 SkyTaper files are used in a conventional torque-limited endodontic contra-angle or endodontic motor. F6 SkyTaper files are supplied sterile-packed and are ready for use.*

Intelligent solution - F6 SkyTaper

- *Root canal preparation in rotary mode with only one file of taper .06.*
- *Complete range of files in 5 sizes to assure adaptation to all root canal anatomies.*
- *Thanks to their well thought-out design, the files are extremely flexible yet efficient.*
- *Sterile-packed single-use files*

Root canal preparation with only one file

- 1 *Ream the canal entrance area to a sufficient size, for example using the Opener ref. OP10L19.*
- 2 *Create a glide path, for example with a K file in size 010.*
The appropriate F6 SkyTaper file is chosen depending on the size of the last (manual) file used. Skip one ISO file size. Skip one ISO size.
- 3 *Root canal preparation with the F6 SkyTaper file*

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
F6 SkyTaper file, taper 06, ring-shaped colour code

Single-use file, sterile packed, for canal preparation to the full working length in picking motion, nickel-titanium

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F360® - Anything other than complicated

- Root canal preparation with only 2 files
- Each file is used in rotary mode to its full working length
- Single-use files
- Sterile-packed

Applied intelligence

Prior to using the F360 files, an Opener (e.g. AK10L19) is used to ream the canal access zone to a larger size and remove all infected substance. Subsequently, patency of the canal has to be ensured by means of suitable manual instruments. In most cases, the root canal can now be prepared with just 2 files. The red F360 file in size 025 is first in line, followed by the green F360 file in size 035. The use of the first F360 file in size 025 may be sufficient in case of very narrow canals. On the other hand, the use of additional files in sizes 045 (white) and 055 (red) may be required in case of wide canals.

The F360 files are used in rotary "picking" motion, i. e. they are moved along the entire working length in a dabbing, picking manner so that the files are kept moving at all times.

Thanks to their taper of 04, all F360 files are extremely flexible so that they are perfectly able to adapt to different root canal morphologies.

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

EndoPilot mobil

The EndoPilot mobil is an endodontic motor and apex locator all in one, ensuring an efficient and safe preparation of the root canal.

This torque and speed controlled endodontic motor is provided with coloured LED lights indicating the direction of rotation, the torque limit or the position of the apex. These features ensure an efficient preparation of the root canal. Thanks to the fully insulated electric contra-angle, the integrated apex locator allows an exact, real-time determination of the length because the operator is always in full control of the actual position of the file.

The EndoPilot mobil is provided with a file library containing the characteristics of the F6 SkyTaper and F360® files and other commonly used file systems. The "MyFile" function enables the operator to develop individually adapted sequences and to save them.

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The EndoPilot mobil is provided with a large, easy to read touch screen which guides the operator through the menu. All functions are found quickly and easily.

Thanks to the reduced size of the retainer and the wireless, radio-controlled foot switch, the battery-operated EndoPilot mobil is particularly easy to handle and adapts perfectly to the requirements in your dental practice.

Advantages:

- Endodontic motor and apex locator all in one
- Fully insulated electric motor and contra-angle
- File library with many preset file systems
- Option to develop individual sequences (MyFile)
- Torque limiting function, equally effective in reverse operation
- Apex locator with impulse measurement
- The preparation length is manually variable at the apex locator
- Once the apex has been reached, the motor stops instead of switching to reverse operation, which would transport the chips towards the apex
- Torque reduction in the vicinity of the apex
- Battery-operated
- Wireless, radio-controlled foot switch
- Clearly laid out menu with easy-to-use touch screen
- Can be updated to allow for future technical developments

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EP0013 EndoPilot mobil

E-Drive - Torque limited endodontic contra-angle

The E-Drive contra-angle is directly placed on the coupling of the micro motor. It can be used with all commonly used nickel titanium file systems (e.g. F6 SkyTaper and F360®). Due to the fact that the torque can be transmitted at 5 different levels and the E-Drive can be directly connected to an apex locator, the mechanical preparation of the root canal can take place in complete safety. Should the file jam in the canal because the torque has been exceeded, it can be released by retro rotation (left-right movements).

Advantages:

- The E-Drive can be placed directly onto the coupling of the micro motor
- The torque can be transmitted at 5 different levels
- Optional setting of intermediate levels
- The torque can be set directly at the selection ring of the contra-angle
- The E-Drive can be directly connected to an apex locator
- The apex can be located without using a file clamp
- Transmission 115:1
- Maximum speed of the motor: ω_{max} 40,000 rpm
- Performs full rotations 360°
- Safety thanks to retro rotation
- Small head to ensure unobstructed view
- Can be sterilised in the autoclave at up to 135°C

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



EasySeal - Root filling material on the basis of epoxy resin

After a successful treatment of the root canal, the root filling is responsible for tightly sealing the canal to prevent reinfection.

A permanent apical seal of the root canal can be achieved with EasySeal, a root filling material based on epoxy resin. EasySeal is particularly suitable for orthograde filling techniques. It is selfsterile to prevent reinfection, radiopaque and dimensionally stable. EasySeal comes in a double-chamber syringe for easy, safe handling without need of previous mixing. One syringe contains 12 g which is enough for approx. 20 root fillings. The workable time of EasySeal is 15 minutes. Full setting is reached at 37°C body temperature within 24 hours.

Root filling materials based on epoxy resin, such as EasySeal, are classed as gold standard thanks to their outstanding properties, for example excellent adhesion to dentin, low solubility and high thermo stability.

Leakage test with dye *

- 1. EasySeal (no dye penetration)**
- 2. Commercially available sealer X (no dye penetration)**
- 3. Commercially available sealer Y (dye penetration)**

** Natural tooth, prepared up to ISO 035, filled with different sealers and gutta-percha, immersed in dye penetrant for 48 hours, removed from the dye penetrant, cut in thin slices in intervals of 2-6 mm from the apex and examined for residues of the dye penetrant under a microscope.*

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F360 Fill

F360 Fill

F360 Fill is a carrier-based filling system for thermoplastic, three-dimensional root canal fillings. The F360 Fill obturators consist of a plastic core coated with thermoplastic gutta-percha. The obturator is heated in the F360 Fill oven in order to guarantee a tight obturation of the root canal. F360 Fill is a universal post system that ideally complements for example the F360® files.

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
Guttapercha cutter

For cutting guttapercha or the plastic carrier in case of using a carrier-based obturation system

To be used preferably in the red contra-angle without cooling agent, applying low contact pressure

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Endo Rescue Kit - Removal of instrument fragments

The fracture of an instrument during an endodontic treatment not only causes the dentist enormous stress, it also poses an increased risk of post-endodontic complications to the patient.

The removal of such fragments is often extremely difficult and almost impossible to plan.

The Endo Rescue Kit provides a simple and systematic solution, providing access to the opening of the root canal and allowing the removal of the fractured instrument. Once straight access to the fragment has been created by means of a conventional endodontic drill and two Gates burs, two specifically developed instruments greatly simplify a previously complicated procedure:

A centre drill (RKP) exposes the coronal part of the fragment. An extremely fine trepan bur (RKT) is then placed onto the fragment which is seized by the bur and held in place by dentin residues. The fragment is then pulled out of the root in an anticlockwise direction. Thanks to the small diameter of these two instruments, the fragment can be extracted whilst removing the absolute minimum of dentin.

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
150.155.000


Clinical sequence

1. Initial situation. Presence of a fragment of a fractured instrument in a mesial root.
2. Recreation of the access cavity (H269GK.315.016).
3. The coronal curvature is straightened with axial movements to create direct access to the fragment (G180A.204.110).
4. Preparation of the access right down to the fragment (G180A.204.090).
5. Exposure of the fractured surface by drilling around the fragment (RKP.204.090).
6. The fragment is surrounded and seized. Anti-clockwise rotation. Removal of the fragment which is firmly held in the trepan bur by the residues of dentin (RKT.204.090).
7. The last third can now be correctly prepared.
8. Final situation after filling.

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

Endo Rescue Trepan bur

To be used in anticlockwise rotation

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

Endo Rescue Kit

For the removal of fractured instruments

					
	H269GK.315.016	1			
	G180A.204.110	1			
	G180.204.090	1			
	RKP.204.090	1			
	RKT.204.090	1			
	155.000.	1			

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

Gutta-percha remover without cutting edges

Plastification of gutta-percha due to frictional heat caused by rotation, nickel-titanium alloy

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


Root filler "L", stainless spring steel

In countries other than Germany and Austria the packing unit is 4 instead of 6

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Trepan bur for exposing fragments in the root canal

9114

	12
--	----

9114.000. ...

Stop per strumenti canalari, silicone
Stopper for root canal preparation instruments, silicone

9138

	1	1	1	1	1	1	1	1
Misura - Size	1	2	3	4	5	6	7	8

9138.000. ...

1 2 3 4 5 6 7 8

Dispenser Endo Stop con 100 stop in silicone
Endo Stop Dispenser with 100 silicone stoppers

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MULTIflex® is a registered trademark of KaVo Dental GmbH, Germany
Komet Italia S.r.l.