Endodonzia *Endodontics*

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Brasseler®, Komet®, Art2®, CeraBur®, CeraCut®, CeraDrill®, CeraFusion®, CeraPost®, DC1®, DCTherm®, Derminator®, FastFile®, F360®, F6 SkyTaper®, H4MC®, OccluShaper®, OptiPost®, PolyBur®, PrepMarker®, Procodile®, R6 ReziFlow®, TissueMaster®, TMC®, TissueMaster Concept® e Visio-Soft® sono marchi registrati della Gebr. Brasseler GmbH & Co. KG.

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Sirona® is a registered trademark of the company Dentsply Sirona. PIEZOlux™, MULTIflex™ and SONICflex™ are trademarks of the company KaVo. Proxeo®, Synea® and Alegra® are registered trademarks of the company W&H.

Brasseler®, Komet®, Art2®, CeraBur®, CeraCut®, CeraDrill®, CeraFusion®, CeraPost®, DC1®, DCTherm®, Derminator®, FastFile®, F360®, F6 SkyTaper®, H4MC®, OccluShaper®, OptiPost®, PolyBur®, PrepMarker®, Procodile®, R6 ReziFlow®, TissueMaster®, TMC®, TissueMaster Concept® and Visio-Soft® are registered trademarks of Gebr. Brasseler GmbH & Co. KG.

Some of the products and designations mentioned in the text are trademarked, patented or copyrighted. The absence of a special reference or the sign ® should not be interpreted as the absence of legal protection.

This publication is copyrighted. All rights, also with regard to translation, reprint and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or reprocessed using electronic systems in any form or by any means (photocopying, microfilm or other methods) without the written permission of the editor.

Colours and products are subject to alterations. Printing errors excepted.

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78 Portastrumenti
Insert tray

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Access cavity

The creation of a perfect access cavity is the first step of any endodontic treatment. In many cases, this is more complicated than the subsequent preparation of the root canal. The formation of the access cavity is divided into two stages: the preparation of the primary cavity (access to pulp chamber) and the secondary cavity (access to the root canal system). The success of the endodontic treatment largely depends on these two vital first steps - the creation of sufficient space and unobstructed vision. Komet offers a comprehensive range of special instruments for this.

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Endo reamer with safe end

EndoGuard The course is set for Endo success

The creation of a perfect access cavity is the foundation of any successful endodontic treatment.

Used immediately after opening the pulp chamber, the EndoGuard helps to carry out this important treatment step with optimum efficiency and in perfect safety.

- Removal of dentin overhangs for improved vision and easier detection of the canal entrances
- The straight access to the root canal minimizes the risk of inadvertent canal transportation and fracture of the file
- Smooth, non-cutting tip to protect the floor of the pulp chamber
- Tapered shape to prevent the preparation of undercuts
- Efficient cross-cut toothing for perfectly controlled, yet smooth work with minimum vibration

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Combination bur with diamond coated working part and carbide tip, for lateral substance removal when preparing an access cavity

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For creation of an endo access cavity

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
EndoTracer

The EndoTracer is an endodontic instrument specially intended for the preparation of the endodontic access cavity, notably for the preparation of isthmuses.


- Thanks to their special design, the petite round burs - especially the sizes 004 and 006 - are ideally suitable for the fine shaping of isthmuses and canal entrances
- The EndoTracer comes in two lengths and in 6 different sizes. This means that the range comprises the perfect instrument for any clinical situation
- Provided with a particularly long, slim neck, the EndoTracer allows unobstructed view past the instrument into the access cavity
- With its total length of 34 mm and an extra 3 mm in neck region, the EndoTracer is now even more suitable for work under a microscope

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EndoTracer for the preparation of the endodontic access cavity,
notably for the preparation of isthmuses

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				
	H1SML31.205.006	1		
	H1SML31.205.008	1		
	H1SML31.205.010	1		
	H1SML31.205.012	1		
	H1SML31.205.014	1		
	H1SML34.205.006	1		
	H1SML34.205.008	1		
	H1SML34.205.010	1		
	H1SML34.205.012	1		
	H1SML34.205.014	1		

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
EndoExplorer **Form follows function**

Developed in close cooperation with the endodontic specialist Dr. Hans-Willi Herrmann, this new instrument set is designed for the ergonomic, tissue-friendly primary and secondary preparation of the endodontic access cavity.

- The instrument design is optimally adapted to the requirements of microscope users
- Petite instrument head in combination with a long, slender neck for unobstructed view
- Controlled, precise work almost without pressure thanks to the sharp toothing
- Tapered instrument head to allow controlled guidance of the instrument and almost non-invasive work
- Completely made of tungsten carbide for maximum runout accuracy even after several uses

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

EndoExplorer for the preparation of the endodontic access cavity, tungsten carbide

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
Pulp bur "M ller", stainless steel

Coronal reaming

The opening of the root canal system is followed by the preparation of the root canals. Given that most bacteria are found in the coronal third of the root canal, it is particularly important to clean this section of the canal right at the start, making sure that the subsequently used files do not spread the bacteria towards the apex. This step has the added advantage of facilitating the work of the next files and improving the dentist's field of vision. Komet provides a vast line of instruments for this, from the well-established Gates-Glidden burs to openers. These files for the mechanical preparation of root canals are provided with a short working part and large taper to facilitate particularly efficient work.

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STERILE R

Opener

This Opener is ideally suited for the fast and generous enlargement of the canal entrance zone. It reliably removes most bacteria from the canal right at the start of the treatment. The fast-cutting design of the file ensures quick and thorough removal of the infected coronal tissue. What's more, the construction of the opener opens new perspectives: an improved view of the work field, especially in combination with magnifying glasses or a microscope. The preparatory work done by the Opener relieves the strain on the subsequently used file system and creates ideal conditions for the preparation further down the root canal.

By the way, it does not matter if your prefer working to a reciprocating or rotary method: both versions of the Opener are suitable for universal use and can be combined with any file system and preparation technique.

- Quick enlargement of the canal entrance zone
- Suitable for universal use
- Fast-cutting file design
- Less strain on the subsequently used file systems
- Improved vision
- Sterile

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Opener, taper .08, length 19 mm, sterile packed, for the straight root canal entrance area with small lumen, nickel-titanium

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
Opener, taper .10, length 15 or 19 mm, sterile packed, for the straight root canal entrance area, nickel-titanium

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*Reamer for the root canal entrance area, taper .10, length 15 or 19 mm, nickel-titanium with TiN coating
In countries other than Germany and Austria the packing unit is 5 instead of 6*

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
Reamer Gates Glidden "G", stainless steel

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


Creation of a glide path

The creation of a glide path is indispensable for ensuring successful endodontic treatments. This step facilitates the correct estimation of the anatomy of the root canal, thus ensuring that the subsequently used files for root canal preparation work safely and efficiently. The glide path can be created manually or mechanically – Komet offers a complete solution for both options.

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PathGlider

PathGlider

The PathGlider is intended for the mechanical creation of a glide path. Made of highly flexible nickel titanium and provided with a non-cutting instrument tip, the PathGlider perfectly follows the course of the canal, thereby not only reducing the risk of inadvertent canal transportation and formation of undesirable steps, but also saving precious time. Thanks to its taper of .03, the instrument clears the way for the subsequently used mechanical files (e.g. F360: taper .04), making sure that they work safely and efficiently. The PathGlider can be combined with any file system. Its unique kite-shaped cross-section guarantees smooth canal walls and excellent control of the file inside the canal.

Advantages over an exclusively manual creation of a glide path:


- Reduced risk of inadvertent transportation of the canal and creation of undesirable steps
- Economy of time
- Unique kite-shaped cross-section for smooth canal walls and excellent control of the file inside the canal
- Transition taper .03 for a gradual, gentle transition to the subsequently used mechanical file systems
- Sterile

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PathGlider file, taper .03, single-use file, sterile packed, for the mechanical establishment of a glide path, nickel-titanium

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
Manual file for probing and checking patency of root canals, stainless, heat-strengthened steel

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Manual preparation

Three basic types of files are available for the manual preparation of the root canal: K-burs, K-files and Hedstroem files. The preparation of root canals is governed by the standard ISO 3630 which applies to root canal instruments. In compliance with ISO provisions, all manual files come with a taper of .02 and a toothed working part with a length of 16 mm. Beyond the requirements of the standard, manual files made by Komet offer the following additional advantages:

- Radiopaque stoppers preassembled at our works
- Clearly visible depth marks
- Ergonomic handle with knobs to prevent slipping
- Handle provided with a hole to insert a safety chain

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



Instruments for root canal preparation, root canal reamers "K", stainless spring steel

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
Instruments for root canal preparation, "K" files, stainless spring steel

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

Instruments for root canal preparation, Hedstroem files „H“, stainless spring steel

Sortimente:
Assortments:

L = 21 mm


174 21.654.S2

1 x 045 - 080

174 25.654.S2

174 31.654.S2

17421
17425
17431

		6	6	6	6	6	6	6
Misura - Size	Ø 1/100 mm	045	050	055	060	070	080	090
Impugnatura - Handle								
340 654 650453 ...								
17421.654. ...		045	050	055	060	070	080	-
340 654 651453 ...								
17425.654. ...		045	050	055	060	070	080	090
340 654 653453 ...								
17431.654. ...		045	050	055	060	070	080	-

Lime di Heström tipo «H», acciaio inossidabile elastico
Instruments for root canal preparation, Hedstroem files "H", stainless spring steel


Mechanical preparation

Apart from the basic task of removing as much infected tissue from the root canals as possible, the mechanical preparation of root canals is focused on providing the canal with a tapered shape. This allows reliable disinfection and tight obturation of the canal. With NiTi instruments, high-quality work can be done in a reliable manner. In the last few years, several power systems for mechanical files have established on the market. Komet's comprehensive range comprises rotary file systems and a system for reciprocating use.

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F6 SkyTaper

A new level of flexibility

F6 SkyTaper, the single file system with taper .06 files made of nickel-titanium which offers uncompromising flexibility to both endo specialists and general dentists.

- Rotary preparation along the full working length with one taper .06 file
- All-comprehensive range of files: NiTi files in five sizes and three lengths to suit all root canal anatomies
- Efficient double-S cross-section for thorough cleaning
- Excellent preservation of the course of the canal
- Unrestrictive choice of the obturation method thanks to taper .06
- Rotary use in all torque-limited endodontic contra-angles and motors
- Sterile packed single-use files

Step by Step

① Generous preparation of the entrance of the root canal, e.g. with the opener OP10L19.

② Creation of a glide path, for example with a K-file in size 015.

The suitable F6 SkyTaper file is chosen according to the size of the previously used (manual) file.

Rule: Skip one ISO file size.

③ Preparation of the root canal with the F6 SkyTaper file.

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
F6 SkyTaper file, taper .06, single-use file, sterile packed, for canal preparation to the full working length in pecking motion, nickel-titanium

28

new

GPF 06

		100	100	100	100	100
		Yellow	Red	Blue	Green	Black
Misura - Size	Ø 1/100 mm	020	025	030	035	040

GPF06.000. ...	020	025	030	035	040
----------------	-----	-----	-----	-----	-----

F6 SkyTaper punte di guttapercha

Con codifica colori, graduate e radioopache

Lunghezza 28 mm

F6 SkyTaper Guttapercha points

Colour coded, graduated and radiopaque

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
F360 - Anything other than complicated

- Root canal preparation with only 2 files
- Each file is used in rotary mode to its full working length
- Single-use files
- Sterile-packed

Applied intelligence

Prior to using the F360 files, an Opener (e.g. OP10L19) is used to ream the canal access zone to a larger size and remove all infected substance. Subsequently, patency of the canal has to be ensured by means of suitable manual instruments. In most cases, the root canal can now be prepared with just 2 files. The red F360 file in size 025 is first in line, followed by the green F360 file in size 035. The use of the first F360 file in size 025 may be sufficient in case of very narrow canals. On the other hand, the use of additional files in sizes 045 (white) and 055 (red) may be required in case of wide canals.

The F360 files are used in rotary "pecking" motion, i. e. they are moved along the entire working length in a dabbing, pecking manner so that the files are kept moving at all times.

Thanks to their taper of .04, all F360 files are extremely flexible so that they are perfectly able to adapt to different root canal morphologies.

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

Well thought-out in every detail

- Manageable file sequence: only 2 files are needed to suit most canal anatomies
- Use in rotary mode in a conventional torque-limited endodontic contra-angle or motor
- All files are used with an identical torque of 1.8 Ncm so that there is no need to adjust the torque of the endodontic contra-angle when changing the file
- Simple use according to the single length principle
- Optimum respect of the original canal anatomy
- Excellent cleaning result and quick mechanical preparation thanks to the innovative instrument design: cross-section in the shape of a double S for increased flexibility, large chip spaces to flush out all debris, blade twist with dynamic increase
- Reduced risk of instrument fracture and no risk of cross-contamination as the F360 files are intended for single use only
- No need to clean, disinfect and sterilise instruments. No need to document and control the number of uses as the F360 single-use files are sterile packed
- Flexible files made of NiTi with small taper .04 assure optimal adaptation to nearly all canal anatomies
- The files are also available in the additional sizes 045 and 055 and come in 3 lengths (L21, L25 and L31) to allow for optimal preparation of all canal anatomies

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F360 file, taper .04, single-use file, sterile packed, for canal preparation to the full working length in pecking motion, nickel-titanium

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
Procodile

Procodile

- Unrivalled flexibility even in curved root canals thanks to an innovative, variably tapered file core
- Greater efficiency during preparation thanks to an enlarged chip space
- Increased patient safety thanks to greater resistance to cyclic fatigue
- Suitable for use in all common left-cutting reciprocating power units and with the patented ReFlex movement in the EndoPilot

Step by Step

- 1 Generous opening of the canal entrance area for example with the reciprocating Opener OPR08L19
- 2 Creation of a glide path, for example with a K-file in size 015. The appropriate Procodile file is chosen depending on the size of the previously used (manual) file
General rule: Skip one ISO size
- 3 Preparation of the root canal with the Procodile file



Obturation with a matching guttapercha point

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

R6 ReziFlow Reciprocating treatments have never been so gentle

The design of the R6 ReziFlow is optimally adapted to reciprocating movements. Provided with a reduced number of spirals, the file has a low screw-in effect which allows gentle, yet effective treatments. The operator stays in complete control of the file throughout the session. At the same time, the intelligent cross-section of the file ensures effective substance removal whilst perfectly preserving the original curvature of the canal.

- One file reciprocating system made of NiTi
- Gentle and controlled work, no screw-in effect
- Fast cutting, perfectly flexible double-S cross-section
- Cuts to the left, therefore suited for use in all established reciprocating motors

Step by Step

Generous opening of the canal entrance area for example with the Opener OP10L19.

- 1 Creation of a glide path, for example with a K-file in size 015. The appropriate R6 ReziFlow file is chosen depending on the size of the previously used (manual) file.
General rule: Skip one ISO size.
- 2 Reciprocating preparation of the root canal with the R6 ReziFlow file.
- 3 Obturation with a matching gutta-percha point.

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
R6 ReziFlow, size 025-040 with taper .06, size 045-050 with taper .04, single-use file, sterile packed, for the reciprocating preparation of root canals in a "pecking motion" along the entire working length, nickel titanium

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Rinsing and activation

About a third of the root canal surface is not reached by mechanical cleaning with files. What's more, mechanical cleaning is not sufficient to reduce the bacterial count enough to create good healing conditions. The chemical disinfection of the endodontium is therefore essential for the successful preparation of the root canal. The activation of the rinsing solution can notably increase the effectiveness of chemical cleaning. Komet offers a NiTi tip specially designed for sound activated rinsing of the root canal.

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SonicLine

Sonic tip for activating rinsing liquids during endodontic treatments. Thanks to sound activated movements and micro currents, the efficiency of rinsing liquids inside the root canal is greatly increased, safely removing bacteria, residues of pulp tissue, loose dentin chips and the smear layer.

Advantages:

- Even more efficient rinsing of the root canal
- Made of highly flexible nickel titanium, with a titanium-nitride coating
- Small instrument diameter for rinsing narrow root canals
- The same instrument diameter for all canals, thus eliminating the need to change the instrument
- Safe end and absence of toothing to avoid inadvertent removal of substance from the canal wall
- Laser marks to indicate the depth

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
Obturation

After the preparation and chemical cleaning of the root canal system, the canals have to be tightly sealed with a suitable obturation to prevent reinfection. Komet provides a comprehensive range of paper and gutta-percha points as well as F360 Fill, a carrier-based filling system. The line is completed by Easy Seal, a sealer based on epoxy resin.

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EasySeal - Root filling material on the basis of epoxy resin

After a successful treatment of the root canal, the root filling is responsible for tightly sealing the canal to prevent reinfection.

A permanent apical seal of the root canal can be achieved with EasySeal, a root filling material based on epoxy resin. EasySeal is particularly suitable for orthograde filling techniques, it is radiopaque and dimensionally stable. EasySeal comes in a double-chamber syringe for easy, safe handling without need for previous mixing. The workable time of EasySeal is 15 minutes. Full setting is reached at 37°C body temperature within 24 hours.

Root filling materials based on epoxy resin, such as EasySeal, are classed as gold standard thanks to their outstanding properties, for example excellent adhesion to dentin, low solubility and high thermo stability.

Leakage test with dye*

1. EasySeal (no dye penetration)
2. Commercially available sealer X (no dye penetration)
3. Commercially available sealer Y (dye penetration)

** Natural tooth, prepared up to ISO 035, filled with different sealers and gutta-percha, immersed in dye penetrant for 48 hours, removed from the dye penetrant, cut in thin slices in intervals of 2–6 mm from the apex and examined for residues of the dye penetrant under a microscope.*

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EasySeal

Permanent root filling material on the basis of epoxy polymer

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20 punte Endo EasySeal
20 EasySeal endo tips

F360 Fill

F360 Fill is a carrier-based filling system for thermoplastic, three-dimensional root canal fillings. The F360 Fill obturators consist of a plastic core coated with thermoplastic gutta-percha. The obturator is heated in the F360 Fill oven in order to guarantee a tight obturation of the root canal. F360 is a universal carrier based post system that ideally complements our file systems F360, F6 SkyTaper and R6 ReziFlow.

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Clinical sequence

1. Prepared, rinsed and dried root canal
2. Determination of the size with the F360 Fill Verifier
3. Choose suitable F360 Fill Obturator and heat up in the oven (green light and + acoustic signal when heated)
4. Apply sealer (e.g. EasySeal 9978)
 - with F360 Fill Verifier or paper point
 - Apply a very thin coat of sealer onto the canal walls
 - Unsuitable for Lentulo
5. Insert obturator into the root canal up to the working length, applying low pressure
6. Remove the protruding core of the obturator and any excess gutta-percha
 - Allow gutta-percha to cool for approx. 3-4 minutes
 - Break off handle with plastic core by hand
 - or with a gutta-percha cutter

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
Guttapercha cutter

For cutting guttapercha or the plastic carrier in case of using a carrier-based obturation system

To be used preferably in the red contra-angle without cooling agent, applying low contact pressure

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
Root filler "L", stainless spring steel

In countries other than Germany and Austria the packing unit is 4 instead of 6


Retreatment

In case of a failed endodontic treatment, for example caused by persistent bacteria or a recolonization of the endodontium, a retreatment is often the only chance to save the tooth. The objectives of an orthograde retreatment are to remove the existing filling material to the extent possible, to access parts of the canal that had not yet been reached and to guarantee a complete re-disinfection of the root canal. The range provided by Komet to this end includes instruments for removing all types of different filling material and a special instrument set for removing fractured files that are stuck inside the root canal.

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Torque: 1,8 Ncm

STERILE

Endo ReStart A new era in retreatment

The goal of every retreatment is to correct an existing endodontic therapy. Thanks to the NiTi retreatment files contained in the Endo ReStart system, this task can now be tackled effectively and under perfect control in just a few steps.

- Endo ReStart is a new NiTi file system for efficient retreatments with just 1-2 instruments
- Thanks to their special blades with dynamic twist, the files of the Endo Restart system ensure effective removal of root fillings
- Thanks to its „safe activity“ instrument tip, we can offer an optimized file to users. The new tip is provided with a tapered, non-cutting end permitting easy penetration of the filling material
- The constant taper of .05 ensures a perfect balance between stability and flexibility

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


Removal of Gutta-percha with
the GPR

- 1. Insert the gutta-percha remover GPR between the gutta-percha and the canal wall, rotating the instrument at an optimal speed of 2000 rpm (without irrigation)
- 2. Once the gutta-percha has softened, increase the speed to the maximum speed of 4000 rpm and insert the instrument further towards the apex
- 3. Pull the gutta-percha remover out of the canal. The plasticized gutta-percha adheres to the instrument

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

Removal of F360 Fill Obturators

Expose the carrier/plastic core at the coronal part of the tooth, i.e. remove the gutta-percha at the top, for example with an EndoTracer in size 010 at an optimum speed of $\odot_{opt}$ 5000 rpm. After that, the carrier can be pulled out with small pincers or special tweezers for removing fragments (ref. 215 or 216).

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EndoTracer for the preparation of the endodontic access cavity, notably
for the preparation of isthmuses

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Endo Rescue Kit - Removal of instrument fragments

The fracture of an instrument during an endodontic treatment not only causes the dentist enormous stress, it also poses an increased risk of post-endodontic complications to the patient. The removal of such fragments is often extremely difficult and almost impossible to plan.

The Endo Rescue Kit provides a simple and systematic solution, providing access to the opening of the root canal and allowing the removal of the fractured instrument. Once straight access to the fragment has been created by means of a conventional endodontic drill and two Gates burs, two specifically developed instruments greatly simplify a previously complicated procedure:

A centre drill (RKP) exposes the coronal part of the fragment. An extremely fine trepan bur (RKT) is then placed onto the fragment which is seized by the bur and held in place by dentin residues. The fragment is then pulled out of the root in an anticlockwise direction. Thanks to the small diameter of these two instruments, the fragment can be extracted whilst removing the absolute minimum of dentin.

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1

Clinical sequence

1. Initial situation. Presence of a fragment of a fractured instrument in a mesial root.

2. Recreation of the access cavity (H269GK.315.016).

3. The coronal curvature is straightened with axial movements to create direct access to the fragment (G180A.204.110).

4. Preparation of the access right down to the fragment (G180A.204.090).

5. Exposure of the fractured surface by drilling around the fragment (RKP.204.090).

6. The fragment is surrounded and seized. Anti-clockwise rotation. Removal of the fragment which is firmly held in the trepan bur by the residues of dentin (RKT.204.090).

7. The last third can now be correctly prepared.

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


Post Endo

Root posts are intended for the reconstruction of teeth after successfully concluded endodontic treatments. It is their task to ensure permanent retention of the coronal build-up. What's more, the weakening of the root during preparation and load transmission has to be kept to the absolute minimum. The comprehensive range offered by Komet includes for example the ER system which has been a mile stone in the post endodontic treatment sector for almost 30 years. It comprises tapered, passive root posts that are available in various different material perfectly adapted to each individual indication.

New additions to the ER system include the two short posts DentinPost X Coated and TitanPost X Coated made of glass fiber reinforces composite and pure titanium. Thanks to their short shanks of merely 6 mm, the root is hardly weakened at all.

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



Sonic tips

The sonic tips from our SonicLine facilitate endodontic treatments with a gentle, conservative touch. Komet not only provides sonic tips for the orthograde preparation of the pulp chamber, but also sonic tips for the retrograde preparation of the root canal as part of an apicoectomy.

Important notes:

The sonic tips made by Komet can not only be used in Komet's sonic hand-piece SF1LM/S, but also

- *in the Scalers made by the co. W&H (Series Proxeo® ZA-55/L/LM/M/LS and Proxeo® ST ZE-55RM/BC, Series Synea® ZA-55/L/LM/M or series Alegra® ST ZE-55RM/BC)*
- *in the SONICflex™ hand-piece provided by the company KaVo (Series 2000N/L/X/LX or series 2003N/L/X/LX)*
- *in the SIROAIR L of the co. Dentsply Sirona*
- *in the sonic hand-pieces Ti-Max, series S970L/KL/SL, made by NSK*

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
SonicLine

Sonic tips for orthograde preparation of the pulp chamber and preparation of the cervical third of the root canal as part of an endodontic treatment.

Advantages:

- Quick preparation and removal of old root fillings
- Easier retrieval of root canals
- Enlargement of obliterated canals
- Preparation of straight canal access cavities without weakening the crown
- Useful for removing hard root fillings, cements or posts
- Controlled, gentle preparation without steps and protruding material
- Excellent vision

The sonic tips can be used in the air driven sonic hand piece SF1LM/S (air scaler) and reprocessed with the help of a rinse adapter in a Miele washer/ disinfecter.

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Tapered

Retrieval of fine and calcified canals, opening of the upper canal portions during revision

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Tapered

Alternative to the SF67 with a more pronounced angle

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Small bud

Finishing of the access cavity, minimal refreshing of the dentin layer and removal of residues of root canal fillings

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



Tapered

Enlargement of long and wide canals, unblocking of fractured instruments, removal of root fillings made of gutta-percha and soft cements

SonicLine

Sonic tips for retrograde endodontic treatments as part of an apicectomy.

Clinical sequence:

Preparative surgical steps

1. Create the access through the jaw bone and work on the bone in the periapical region. Prepare a small bone window with a repositionable lid to ensure a proper preparation of the canal with the endodontic tips for retrograde treatments.
2. Remove the root apex in a right angle to the tooth axis, for example with the tungsten carbide bone cutter Komet H254E
3. Remove any diseased tissue by means of a sharp spoon curette and stop the resulting bleeding

Retrograde root preparation with SonicLine sonic tips

4. Retrieval of the root canal and exposure of the canal entrance with the sonic tips SF56 (curved to the left) or SF57 (curved to the right)
5. Preparation of the cavity up to a depth of 3 mm with the SF16 (curved to the left) or SF17 (curved to the right)
In case of extremely curved root canals, it might be helpful to open and enlarge the canal with the specially shaped sonic tip SF55
6. Preparation of an undercut in order to prevent loss of the canal filling with the tip SF20 (curved to the left) or SF21 (curved to the right)
7. Finally, fill the root with a material suitable for retrograde fillings

Advantages:

- Minimally invasive treatment without need to prepare a large bone window
- Axial work, even in very crowded conditions
- Easier work thanks to double-angled tips
- Slender tips for excellent vision in all jaw regions
- Simplified preparation of undercuts for permanent retention of the retrograde root filling

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4602.000

Set accessori SF1979 e 566
Set cooling adapter SF1979 for sonic tips and mounting wrench 566

SF1979.000.

1

566.000.

1

Motori Endo

Endo motors

Special Endo motors are indispensable when it comes guaranteeing a safe mechanical preparation of the root canal. These motors have to meet special requirements, such as reduced speeds and torque control. Komet's endodontic range comprises a torque reduced endodontic contra-angle and a complete endodontic motor capable of combining preparation and electronic determination of the length.

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EndoPilot

EndoPilot

The EndoPilot is an endodontic motor and apex locator all in one, ensuring an efficient and safe preparation of the root canal.

This torque and speed controlled endodontic motor is provided with coloured LED lights indicating the direction of rotation, the torque limit or the position of the apex. These features ensure an efficient preparation of the root canal. Thanks to the fully insulated electric contra-angle, the integrated apex locator allows an exact, real-time determination of the length because the operator is always in full control of the actual position of the file.

The EndoPilot is provided with a file library containing the characteristics of the F6 SkyTaper, F360 and Procodile files and other commonly used file systems. What's more, the operator has the option to develop individually adapted sequences.

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EP0014 EndoPilot

The EndoPilot is provided with a large, easy to read touch screen which guides the operator through the menu. All functions are found quickly and easily.

Thanks to the reduced size of the retainer and the wireless, radio-controlled foot switch, the battery-operated EndoPilot is particularly easy to handle and adapts perfectly to the requirements in your dental practice.

Hardware:

- Endodontic motor and apex locator all in one
- Modern and attractive design
- Fully insulated electric motor and contra-angle
- 7 inch colour touch-screen display
- Wireless, radio-controlled foot switch
- Can be updated by means of micro SD card to allow for future developments
- Battery-operated
- High-quality metal holder and concealed cable routing

Software:

- File library with many preset file systems
- Option to develop individual sequences
- Clearly laid out menu with easy-to use touch screen
- Precise real-time determination of the length thanks to the patented impulse measuring technique
- The preparation length is manually variable at the apex locator
- The motor stops as soon as the preparation length is reached
- Torque reduction in the vicinity of the apex
- Patented ReFlex motion for use in combination with Procodile files
- ReFlex smart provides extra safety thanks to its particularly sensitive reaction to torsional loads
- ReFlex dynamic works particularly fast, thus offering even greater efficiency during preparation

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EndoPilot

Endodontic motor and apex locator in one

E-Drive - Torque limited endodontic contra-angle

The E-Drive contra-angle is directly placed on the coupling of the micro motor. It can be used with all commonly used nickel titanium file systems (e.g. F360 and F6 SkyTaper). Due to the fact that the torque can be transmitted at 5 different levels, the mechanical preparation of the root canal can take place in complete safety. Should the file jam in the canal because the torque has been exceeded, it can be released by retro rotation (left-right movements).

Advantages:

- The E-Drive can be placed directly onto the coupling of the micro motor
- The torque can be transmitted at 5 different levels
- Optional setting of intermediate levels
- The torque can be set directly at the selection ring of the contra-angle
- Transmission 115 :1
- Maximum speed of the motor:
 $\odot_{\text{max.}} 40,000 \text{ rpm}$
- Performs full rotations 360°
- Safety thanks to retro rotation
- Small head to ensure unobstructed view
- Can be sterilised in the autoclave at up to 135°C

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Auxiliaries

To support the dentist with helpful tools before, during and after the treatment, Komet offers specially matched accessories and instrument blocks, for example interim blocks for the intermediate storage of endodontic files during the preparation of the root canal or bur blocks suitable for cleaning, disinfection and sterilization of instruments in the autoclave.

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Insert tray

Endodontic instruments can be stored in a clearly arranged treatment tray.

- For 28 endodontic instruments (handpiece and contra-angle)
- PP, 64 x 59 x 19 mm

Sterilization container made of high-performance plastic (PPSU)

- No seals, no maintenance required. More than 2000 sterilization cycles
- Suitable for vacuum steam sterilization
- ePFTE filter. The long-lasting ePFTE filter remains in the sterilization container during reprocessing (manual or mechanical)
- Transparent material - the contents are visible from the outside
- Two or more containers can be laterally connected

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Wash box

The wash box is intended for mechanical cleaning and disinfection of instruments in the thermo disinfectant. Following an initial pre-cleaning step, the instruments to be cleaned are placed into the insert trays that are part of the chosen treatment system. The loaded insert tray can then be put into the wash box. It is then cleaned and disinfected in the thermo disinfectant.

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Bur block made of stainless steel with 16 blue silicone plugs as universal instrument holders, for FG and RA instruments with a maximum length of 33 mm

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Bur block made of stainless steel with 12 blue silicone plugs as universal instrument holders, for FG and RA instruments with a maximum length of 48 mm

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Bur block made of stainless steel with 12 blue silicone plugs as universal instrument holders, for FG and RA instruments with a maximum length of 68 mm

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MULTIflex™ is a trademark of KaVo
Komet Italia S.r.l.